

Human Granzyme B ELISA Kit

Catalog No: 850.790.192

SPECIFICITY : Recognizes both natural and recombinant human Granzyme B

RANGE : 31.25 pg/ml - 1000 pg/ml

SENSITIVITY : < 20 pg/ml


INCUBATION : 3 hours

SAMPLE TYPES : Serum - Plasma - Cell culture supernatant

SAMPLE SIZE : 100 μ l

CROSS REACTION : No cross reactivity with Trail, TNF- α , IFN- γ , Fas, IL-1 α , IL-4, IL-6, IL-8, IL-13

KIT SIZE : 2 x 96 tests



KIT CONTENT : Pre-coated 12 strip plates, biotinylated secondary antibody, standards, controls (when available), buffers, Streptavidin-HRP, TMB, Stop Reagent.

Use only for research purposes. Not for human use.



Cell Sciences, Inc.
480 Neponset Street
Bldg 12A
Canton, MA 02021

Toll Free: 888-769-1246
Phone: 781-828-0610
Fax: 781-828-0542

E-mail: info@cellsciences.com
Web Site: www.cellsciences.com