

D-glutamate Conjugate

Catalog No.	AG024	Quantity:	1 mg
Description:	D-glutamate Conjugate: D-glutamate-G-Pc (Protein Carrier). G = Glutaraldehyde.		
Applications:	Immunohistochemistry, immunocytochemistry.		
Storage & Handling:	Lyophilized antigen are stable at least 1 year. For reconstitution, use sterile ddH ₂ O. This solution is stable at +4°C for 15 days.		
Related Product(s):	Catalog No. PA2046 : Anti-conjugated D-aspartate (rabbit).		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

