

Mouse Anti-Human CD107a Clone B-T47 Azide Free mAb

Catalog No.	CDM059	Quantity:	200 µg
Alternate Names:	Lysosomal-Associated Membrane Protein 1, LAMP1, LAMPA, LGP120		
Description:	Mouse Anti-Human CD107a Clone B-T47 Azide Free monoclonal antibody		
Concentration:	1.0 mg/ml		
Hybridoma:	Myeloma X63/AG.8653 x Balb/c spleen cells		
Specificity:	Recognizes a 100-120 kDa protein		
Host:	Mouse		
Immunogen:	U937 cell line		
Isotype:	IgG1		
Clone:	B-T47		
Formulation:	0.2 µm sterile filtered solution in PBS. Treated to remove endotoxins. Carrier and preservative free.		
Purification:	Ion exchange chromatography		
Applications:	Flow Cytometry The optimal concentration should be determined by the user for each specific application.		
Storage & Stability:	Stable for 1 year at 2-8°C. DO NOT FREEZE.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

