

SELE

Mouse Anti-Human E-Selectin/CD62E Clone B-S3 Azide Free mAb

Catalog No.	CDM105	Quantity:	200 µg
Alternate Names:	CD62E, ELAM, ELAM1, ESEL, LECAM2, E-selectin, ELAM-1, endothelial adhesion molecule 1; endothelial leukocyte adhesion molecule 1; leukocyte endothelial cell adhesion molecule 2; leukocyte-endothelial cell adhesion molecule 2		
Description:	Mouse Monoclonal Antibody against Human CD62E Clone B-S3. E-Selectin is a protein found in cytokine-stimulated endothelial cells and is thought to be responsible for the accumulation of blood leukocytes at sites of inflammation by mediating the adhesion of cells to the vascular lining. Adhesion molecules participate in the interaction between leukocytes and the endothelium and appear to be involved in the pathogenesis of atherosclerosis.		
Concentration:	200 µg / 200 µl		
Gene ID:	6401		
Specificity:	Recognizes the Endothelial Leucocyte Adhesion Molecule-1 (ELAM-1), E-selectin, a 97 kDa protein		
Host:	Mouse		
Immunogen:	Activated human endothelial cells		
Isotype:	IgG1k		
Clone:	B-S3		
Hybridoma:	Myeloma X63/AG.8653 x BALB/c spleen cells		
Formulation:	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free.		
Purification:	Ion exchange chromatography		
Applications:	Suitable for ELISA		
Storage & Stability:	Stable at 2-4°C for 12 months. For longer storage, freeze aliquots. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences[®]
65 Parker Street
Unit 11
Newburyport, MA 01950

Toll Free: 888-769-1246
Phone: 978-572-1070
Fax: 978-992-0298

E-mail: info@cellsciences.com
Website: www.cellsciences.com