

TNFRSF10A

Mouse Anti-Human TRAILR1/DR4 (Clone B-T35) Azide Free mAb

Catalog No.	CDM113	Quantity:	200 µg
Alternate Names:	APO2; CD261; TRAILR1; TRAILR-1; DR4; tumor necrosis factor receptor superfamily, member 10a		
Description:	The monoclonal antibody recognizes human TRAILR-1, a surface protein which, when engaged by its ligand TRAIL, transduces a cell death signal and induces apoptosis.		
Concentration:	1.0 mg/ml		
Gene ID:	8797		
Conjugate:	Unconjugated		
Specificity:	Recognizes both natural and recombinant TRAIL-R1 (DR4), a 43-46 kDa protein		
Host:	Mouse		
Isotype:	IgG1k		
Immunogen:	Recombinant human TRAIL R1/Fc chimera		
Clone:	B-T35		
Formulation:	Liquid in PBS, Sterile-filtered. Carrier and preservative free.		
Purification:	Ion exchange chromatography		
Applications:	ELISA		
Application Notes:	The optimal concentration should be determined by the user for each specific application.		
Storage & Stability:	Stable at 2-4°C for 12 months. For longer storage, freeze aliquots at -20 to -80°C. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences®
65 Parker Street
Unit 11
Newburyport, MA 01950

Toll Free: 888-769-1246
Phone: 978-572-1070
Fax: 978-992-0298

E-mail: info@cellsciences.com
Website: www.cellsciences.com