


FCGR3A

Mouse Anti-Human CD16 Clone B-E16 PE mAb

Catalog No.	CDM154	Quantity:	100 tests
Alternate Names:	CD16, CD16A, FCG3, FCGR3, FCGR3III, FCR-10, FCRIII, FCRIIIA, GFR3		
Description:	Mouse Anti-Human CD16 Clone B-E16 PE mAb		
Immunogen:	Natural killer cells		
Isotype:	Mouse IgG2a		
Clone:	B-E16		
Hybridoma:	Myeloma X63/AG.8653 x BALB/c spleen cells		
Formulation:	Lyophilized from PBS with 5% BSA and 0.1% sodium azide. Precaution: Sodium azide is a poisonous and hazardous substance which should be handled by trained staff only.		
Purification:	Ion exchange chromatography		
Reconstitution:	Use 10 µl to label 10 ⁶ cells or 100 µl of whole blood		
Conjugation:	PE		
Applications:	Suitable for flow cytometry		
Application Notes:	The optimal concentration should be determined by the user for each specific application.		
Storage & Stability:	Stable at +2-8°C for 12 months For longer storage, freeze aliquots		

A typical staining pattern with the B-E16 monoclonal antibody of lymphocytes



NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences®
65 Parker Street
Unit 11
Newburyport, MA 01950

Toll Free: 888-769-1246
Phone: 978-572-1070
Fax: 978-992-0298

E-mail: info@cellsciences.com
Website: www.cellsciences.com