


## CD34

### Mouse Anti-Human CD34 Class III Clone B-F23 Azide Free mAb

<b>Catalog No.</b>	CDM177	<b>Quantity:</b>	200 µg
<b>Alternate Names:</b>	CD34 antigen, hematopoietic progenitor cell antigen CD34		
<b>Description:</b>	Mouse Anti-Human CD34 Class III Clone B-F23 mAb		
<b>Gene ID:</b>	947		
<b>Specificity:</b>	Recognizes the Gp105-120 antigen, a 105-120 kDa protein		
<b>Immunogen:</b>	KG-1a cell line		
<b>Isotype:</b>	Mouse IgG1		
<b>Clone:</b>	B-F23		
<b>Hybridoma:</b>	Myeloma X63/AG.8653 x BALB/c spleen cells		
<b>Formulation:</b>	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free		
<b>Purification:</b>	Ion exchange chromatography		
<b>Applications:</b>	Suitable for flow cytometry Suitable for Immunohistochemistry Does not trigger KG-1a cell aggregation		
<b>Application Notes:</b>	The optimal concentration should be determined by the user for each specific application.		
<b>Storage &amp; Stability:</b>	Stable at +2-8°C for 12 months DO NOT FREEZE		

A typical staining pattern with the B-F23 monoclonal antibody of KG-1a cell line



NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



**Cell Sciences®**  
65 Parker Street  
Unit 11  
Newburyport, MA 01950

Toll Free: 888-769-1246  
Phone: 978-572-1070  
Fax: 978-992-0298

E-mail: [info@cellsciences.com](mailto:info@cellsciences.com)  
Website: [www.cellsciences.com](http://www.cellsciences.com)