


## TNF

### Mouse Anti-Human TNF- $\alpha$ (Clone B-D9) PE mAb

<b>Catalog No.</b>	CDM278	<b>Quantity:</b>	100 tests
<b>Alternate Names:</b>	Tumor necrosis factor; DIF; TNFA; TNFSF2		
<b>Description:</b>	The monoclonal antibody recognizes TNF- $\alpha$ , a multifunctional proinflammatory cytokine that belongs to the TNF superfamily and is mainly secreted by macrophages. TNF- $\alpha$ can bind to and functions through its receptors TNFRSF1A/TNFR1 and TNFRSF1B/ TNFR. It is involved in the regulation of a wide spectrum of biological processes including cell proliferation, differentiation, apoptosis, lipid metabolism, and coagulation.		
<b>Gene ID:</b>	7124		
<b>Conjugate:</b>	PE		
<b>Specificity:</b>	Recognizes both natural and recombinant human TNF- $\alpha$		
<b>Host:</b>	Mouse		
<b>Isotype:</b>	IgG1		
<b>Immunogen:</b>	Recombinant human TNF- $\alpha$		
<b>Clone:</b>	B-D9		
<b>Formulation:</b>	Lyophilized from PBS with 5 % BSA and 0.1 % sodium azide. Precaution: Sodium azide is a poisonous and hazardous substance which should be handled by trained staff only.		
<b>Purification:</b>	Ion exchange chromatography		
<b>Reconstitution:</b>	Reconstitute with 1 ml deionized water.		
<b>Applications:</b>	Intracellular Flow Cytometry		
<b>Application Notes:</b>	Use 10 $\mu$ l of antibody to label $5 \times 10^5$ cells. The optimal concentration should be determined by the user for each specific application.		
<b>Storage &amp; Stability:</b>	Stable at 2-8°C for 6 months after reconstitution. DO NOT FREEZE.		

A typical staining pattern with the B-D9 monoclonal antibody of PMA and ionomycin activated PBL



NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

