

GNLY

Mouse Anti-Human Granulysin Clone B-L38 Azide Free mAb

Catalog No.	CDM313	Quantity:	200 µg
Alternate Names:	T-cell activation protein 519, 519, D2S69E, LAG-2, LAG2, NKG5, TLA519		
Description:	Mouse Monoclonal Antibody against Human Granulysin Clone B-L38. Granulysin is a member of the saposin-like protein (SAPLIP) family and is located in the cytotoxic granules of T cells, which are released upon antigen stimulation. It is present in cytotoxic granules of cytotoxic T lymphocytes and natural killer cells and has antimicrobial activity against <i>M. tuberculosis</i> and other organisms.		
Concentration:	200 µg / 200 µl		
Gene ID:	10578		
Specificity:	Recognizes both natural and recombinant human Granulysin		
Host:	Mouse		
Immunogen:	Recombinant human Granulysin (NKG5 and LAG-2)		
Isotype:	IgG1k		
Clone:	B-L38		
Hybridoma:	Myeloma X63/AG.8653 x BALB/c lymph node cells		
Formulation:	Phosphate-buffered saline. Sterile filtered through 0.22 µm and treated to remove endotoxins. Carrier and preservative free.		
Purification:	Ion exchange chromatography		
Applications:	Recognizes soluble and membrane antigen		
Storage & Stability:	Stable at 2-4°C for 12 months. For longer storage, freeze aliquots. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

