

Mouse Anti-Human IL-35 Clone B-S35 Azide Free mAb

Catalog No.	CDM327	Quantity:	500 µg
Specificity:	Recognizes both natural and recombinant protein (no cross reaction with IL-12 and IL-27)		
Host:	Mouse		
Immunogen:	EBI3 recombinant protein and IL-35/Fc recombinant protein		
Isotype:	IgG1 Kappa light chain		
Clone:	B-S35		
Hybridoma:	Myeloma X63/AG.8653 x Balb/c node cells		
Formulation:	Phosphate-buffered saline. Sterile-filtered through 0.2 µm. Carrier and preservative free.		
Purification:	Ion exchange chromatography		
Applications:	Suitable for ELISA		
Storage & Stability:	Stable at +2-8°C for 12 months. For longer storage, freeze aliquots.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.