


KLRD1


Mouse Anti-Human CD94 (Clone B-D49) Azide Free mAb

Catalog No.	CDM331	Quantity:	200 µg
Alternate Names:	Killer cell lectin-like receptor subfamily D, member 1		
Description:	The mouse monoclonal antibody recognizes human CD94/ KLRD1, an antigen preferentially expressed on Natural Killer (NK) cells. CD94 is classified as a type II membrane protein because it has an external C terminus. NK cells are a distinct lineage of lymphocytes that mediate cytotoxic activity and secrete cytokines upon immune stimulation. Several genes of the C-type lectin superfamily, including members of the NKG2 family, are expressed by NK cells and may be involved in the regulation of NK cell function.		
Gene ID:	3824		
Concentration:	1 mg/ml		
Conjugate:	Unconjugated		
Specificity:	Recognizes human cell surface CD94		
Host:	Mouse		
Isotype:	IgG1κ		
Immunogen:	IL-2 activated NK cells after tolerization with activated PBL (CD56 negative).		
Clone:	B-D49		
Formulation:	Liquid in PBS. Sterile-filtered. Carrier and preservative free.		
Applications:	Flow Cytometry, Immunochemistry		
Application Notes:	IHC: frozen sections The optimal concentration should be determined by the user for each specific application.		
Storage & Stability:	Stable at 2-8°C for up to 1 year. For longer storage, freeze aliquots at -20 to -80°C. Avoid repeated freeze-thaw cycles.		

A typical staining pattern of NK-L cell line with the B-D49 monoclonal antibody.



Labeling with B-D49 (1:200) on frozen section (acetone fixed) of human colon adenocarcinoma. (Background coloration by Hematoxin.)



NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences®
65 Parker Street
Unit 11
Newburyport, MA 01950

Toll Free: 888-769-1246
Phone: 978-572-1070
Fax: 978-992-0298

E-mail: info@cellsciences.com
Website: www.cellsciences.com