


Mouse Anti-Human HLA Class II (DQ) Clone B-K27 Azide Free mAb

Catalog No.	CDM344	Quantity:	200 µg
Description:	Mouse Anti-Human HLA Class II (DQ) Clone B-K27 Azide Free mAb		
Concentration:	200 µg / 200 µl		
Specificity:	Recognizes a 28 kDa and 33 kDa complex		
Host:	Mouse		
Immunogen:	T4+ CTL clone HG-38		
Isotype:	IgG2a		
Clone:	B-K27		
Hybridoma:	Myeloma X63/AG.8653 x BALB/c spleen cells		
Formulation:	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free		
Purification:	Ion exchange chromatography		
Applications:	Suitable for flow cytometry		
Storage & Stability:	Stable at +2-8°C for 12 months For longer storage, freeze aliquots		

A typical staining pattern with the B-K27 monoclonal antibody of lymphocytes



NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

