

## IFNG

### Rat Anti-Mouse Interferon gamma Clone F1 mAb

<b>Catalog No.</b>	CMI032	<b>Quantity:</b>	0.5 mg
<b>Description:</b>	Rat Anti-Mouse Interferon gamma Clone F1 mAb		
<b>Specificity:</b>	Binds with high efficacy to Native and Recombinant Mouse IFN gamma. Poor neutralizing activity.		
<b>Host:</b>	Rat		
<b>Isotype:</b>	IgG <sub>2a</sub>		
<b>Clone:</b>	F3		
<b>Formulation:</b>	Lyophilized in 0.25 mL PBS + 125 mM Trehalose		
<b>Purification:</b>	Ion Exchange Chromatography		
<b>Reconstitution:</b>	<b>Centrifuge vial prior to opening.</b> Dissolve contents by injecting 0.25 mL sterile distilled water.		
<b>Applications:</b>	Coating Antibody in ELISA and ELISPOT assays IHC WB The optimal concentration should be determined by the user for each specific application.		
<b>Storage &amp; Stability:</b>	Store at 2-4°C. 1 year. After reconstitution, store at 2-4°C for one month or for one year at -20°C. Add 0.02% Sodium Azide as a preservative for long term storage. <b>Precaution: Sodium Azide is a poisonous and hazardous substance which should be handled by trained staff only. Avoid repeated freeze-thaw cycles.</b>		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

