

Rat Anti-Mouse IgE Clone 23G3 Biotin mAb

Catalog No.	CMM017	Quantity:	500 µg
Description:	Rat Anti-Mouse IgE Clone 23G3 Biotin mAb		
Concentration:	500 µg/mL		
Specificity:	Mouse IgE		
Immunogen:	Monoclonal anti-DNP mouse IgE		
Isotype:	Rat IgG ₁ κ		
Clone:	23G3		
Conjugate:	Biotin		
Formulation:	Liquid volume of 1.0 mL in a buffer of PBS/NaN ₃		
Isotype Control	Rat IgG ₁ -Biotin		
Applications:	Applications for relevant formats of this clone include ELISA, FLISA, Flow Cytometry, IHC, WB, ELISpot, neutralization.		
Application Notes:	Working dilutions; ELISA BIOT conjugate 1:5,000 - 1: 10,000. Since applications vary, you should determine the optimum working dilution for the product that is appropriate for your specific need.		
Storage & Stability:	Store at 2-8°C. For long term storage at -20°C, aliquot into smaller volumes to avoid multiple freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

