

Mouse Anti-human MMP12 mAb

Catalog No.	CMM110	Quantity:	100 µg
Description:	Mouse Anti-human MMP12 monoclonal antibody		
Alternate Names:	HME, MME, matrix metallopeptidase 12, macrophage metalloelastase, matrix metalloproteinase 12		
Specificity:	Recognizes human MMP12		
Host:	Mouse		
Immunogen:	Recombinant human MMP12		
Isotype:	IgG2		
Formulation:	Lyophilized from a 0.2 µm sterile filtered solution in PBS		
Purification:	Protein G affinity chromatography		
Reconstitution:	Centrifuge vial prior to opening. Add 0.5 ml sterile distilled water.		
Cross-Reactivity:	No cross-reactivity with other human MMPs		
Applications:	Western Blot Immunohistochemistry The optimal concentration should be determined by the user for each specific application.		
Storage & Stability:	Lyophilized antibody is stable for at least 2 years from date of receipt at -20°C to -80°C. Reconstituted antibody can be aliquotted and stored frozen at -20°C to -80°C for at least six months. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

