

## MMP15

### Mouse Anti-Human Matrix Metalloproteinase 15 mAb

<b>Catalog No.</b>	CMM114	<b>Quantity:</b>	100 µg
<b>Alternate Names:</b>	Matrix Metalloproteinase 15, Membrane Type 2 Matrix Metalloproteinase, MT2-MMP, MTMMP2, SMCP-2		
<b>Gene ID:</b>	4324		
<b>Description:</b>	Mouse Anti-human MMP15 monoclonal antibody		
<b>Specificity:</b>	Recognizes human MMP15		
<b>Host:</b>	Mouse		
<b>Immunogen:</b>	Recombinant human MMP15		
<b>Isotype:</b>	IgG1		
<b>Formulation:</b>	Lyophilized from a 0.2 µm sterile filtered solution in PBS		
<b>Purification:</b>	Protein G affinity chromatography		
<b>Reconstitution:</b>	<b>Centrifuge vial prior to opening.</b> Add 0.5 ml sterile distilled water.		
<b>Cross-Reactivity:</b>	No cross-reactivity with human MMP14/16		
<b>Applications:</b>	Western Blot Immunohistochemistry The optimal concentration should be determined by the user for each specific application.		
<b>Storage &amp; Stability:</b>	Lyophilized antibody is stable for at least 2 years from date of receipt at -20°C to -80°C. Reconstituted antibody can be aliquotted and stored frozen at -20°C to -80°C for at least six months. <b>Avoid repeated freeze-thaw cycles.</b>		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

