

Mouse Anti-Rat Kappa Clone K4F5 Biotin mAb

Catalog No.	CMR031	Quantity:	500 µg
Description:	Mouse Anti-Rat Kappa Clone K4F5 Biotin mAb		
Concentration:	500 µg/mL		
Specificity:	Rat Kappa		
Immunogen:	Rat κ hybridoma		
Isotype:	Mouse (BALB/c) IgG _{2a} K		
Clone:	K4F5		
Conjugate:	Biotin		
Formulation:	Liquid volume of 1.0 mL in a buffer of PBS/NaN ₃		
Isotype Control	Mouse IgG _{2a} -Biotin		
Applications:	Applications for relevant formats of this clone include ELISA, FLISA, Immunohistochemistry and Western Blot.		
Application Notes:	Working dilutions; r ELISA 1:5,000 - 1:10,000, for FLISA with FITC conjugate 1:200 - 1:400, for Flow Cytometry ≤ 1 µg/10 ⁶ cells. . Since applications vary, you should determine the optimum working dilution for the product that is appropriate for your specific need.		
Storage & Stability:	Store at 2-8°C. For long term storage at -20°C, aliquot into smaller volumes to avoid multiple freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.