

## Rat Anti-Mouse VEGFR3/FLT4 mAb

<b>Catalog No.</b>	CMV112	<b>Quantity:</b>	100 µg
<b>Alternate Names:</b>	FMS-like tyrosine kinase 4, Chy, Flt-4, VEGFR-3		
<b>Description:</b>	Rat Anti-Mouse VEGFR3/FLT4 monoclonal antibody		
<b>Specificity:</b>	Recognizes mouse VEGFR3		
<b>Host:</b>	Rat		
<b>Immunogen:</b>	Mouse VEGFR3 extracellular domain		
<b>Isotype:</b>	IgG2		
<b>Formulation:</b>	Lyophilized from a 0.2 µm sterile filtered solution in PBS		
<b>Purification:</b>	Protein G affinity chromatography		
<b>Reconstitution:</b>	<b>Centrifuge vial prior to opening.</b> Reconstitute with sterile distilled water.		
<b>Cross-Reactivity:</b>	No cross-reactivity to mouse VEGFR1 and VEGFR2		
<b>Applications</b>	Western Blot Immunohistochemistry  The optimal concentration should be determined by the user for each specific application.		
<b>Storage &amp; Stability:</b>	Store lyophilized antibody at -80°C. Reconstituted antibody is stable for six months in working aliquots at -80°C. <b>Avoid repeated freeze-thaw cycles.</b>		

**NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.**

