


Purified Rabbit Anti-mouse Eotaxin

Catalog No.	CPE001	Quantity:	0.5 mg
Alternate Names:	SCYA11, CCL11, Eotaxin-1		
Description:	For research use only. This polyclonal antibody, which reacts with recombinant and natural mouse eotaxin, was generated using E. coli-expressed amino acid 1-74 of mature mouse eotaxin as an immunogen.		
	<p>Eotaxin is a 74-amino acid eosinophil-chemotactic CC chemokine originally found in bronchoalveolar lavage fluid from allergic inflammatory subjects. It is involved in regulating the recruitment and activation of inflammatory leukocytes, particularly eosinophils. It may play a fundamental role in the development of allergic responses.</p>		
Specificity:	Mouse		
Host:	rabbit		
Immunogen:	Eotaxin		
Isotype:	IgG		
Formulation:	Protein A purified rabbit IgG, 500 µg, with 0.1% sodium azide, lyophilized.		
Reconstitution:	Reconstitute to 1 mg/ml by adding 500 µl PBS.		
Cross-Reactivity:	Reacts with recombinant and natural mouse eotaxin. Cross-reactivity to eotaxin of other species has not been determined.		
Applications:	The tested titer for Western blot is 1:2,000.		
Storage & Stability:	4°C for short term storage or -20°C in small aliquots for long term storage. Avoid repeated freeze and thaw cycles.		



Mouse eotaxin detected in stimulated lung tissue homogenate by Western blot.

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

