

PLG

Native Bovine Plasminogen

Catalog No.	CRP258B CRP258C	Quantity:	1.0 mg 10 mg
Description:	Native Bovine Plasminogen prepared from fresh bovine plasma by immobilized lysine chromatography.		
Concentration:	5.0 mg/ml		
Gene ID:	280897		
Source:	Bovine plasma		
Molecular Weight:	92.0 kDa		
Formulation:	Frozen Liquid in 0.1 M HEPES + 0.1 M NaCl, pH 7.4.		
Purity:	>95% by SDS-PAGE		
Volume:	0.222 ml		
Storage & Stability:	Store in working aliquots at -80°C. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

