

PLAU

Recombinant Human LMW Urokinase FITC

Catalog No.	CRU117B CRU117C	Quantity:	100 µg 1.0 mg
Alternate Names:	Urokinase Plasminogen Activator, PLAU, ATF, UPA, URK, u-PA		
Description:	FITC labeled LMW urokinase purified from auto-catalytic digestion of HMW human urokinase. LMW uPA contains the proteinase domain and is separated from the amino terminal fragment by cleavage at the Lys135-Lys136 bond.		
Concentration:	2.5 mg/ml		
Gene ID:	5328		
Source:	Insect cell culture		
Molecular Weight:	33 kDa		
Formulation:	Frozen liquid in 0.05 M Sodium Acetate + 0.1 M NaCl + 1 mM + EDTA, pH 5.0.		
Purity:	>95% by SDS-PAGE		
Storage & Stability:	Store in working aliquots, protected from light, at -80°C. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

