

Purified Rat IgG2a, K Isotype Control Clone RTK2758

Catalog No.	CSI11905 CSI11906	Quantity:	50 µg 0.5 mg
Description:	Purified Rat IgG2a, K Isotype Control Clone RTK2758		
Purity:	>95% as determined by SDS-PAGE.		
Specificity:	KLH		
Immunogen:	Trinitrophenol + KLH		
Isotype:	Rat IgG2a, K		
Clone:	RTK2758		
Formulation:	0.2 µm filter sterilized liquid in PBS, pH 7.2 with no preservative.		
Purification:	Affinity chromatography		
Endotoxin Level:	<0.1 EU/µg of protein		
Applications:	Flow cytometry		
Application Notes:	The RTK2758 immunoglobulin is useful as an isotype-matched control (for the relevant formats) for Western Blotting, immunoprecipitation, immunohistochemistry, functional assay, immunofluorescence microscopy, immunocytochemistry and immunofluorescent staining (surface or intracellular) for flow cytometric analysis. The optimal concentration should be determined by the user for each specific application.		
Storage & Stability:	Store undiluted at 2-4°C.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.