

## Rat Immunoglobulin G2b kappa Isotype Control Clone RTK4530 Biotin Antibody

<b>Catalog No.</b>	CSI11909 CSI11910	<b>Quantity:</b>	50 µg 200 µg
<b>Description:</b>	The isotype of RTK4530 immunoglobulin is rat IgG2b, κ. This antibody was chosen as an isotype control after screening on a variety of resting, activated, live, and fixed mouse, rat and human tissues.		
<b>Concentration:</b>	0.5 mg/ml		
<b>Host:</b>	Rat		
<b>Immunogen:</b>	Trinitrophenol + KLH		
<b>Isotype:</b>	Rat IgG2b, κ		
<b>Clone:</b>	RTK4530		
<b>Formulation:</b>	Phosphate-buffered solution, pH 7.2, containing 0.09% sodium azide. <b>Precaution:</b> Sodium azide is a poisonous and hazardous substance which should be handled by trained staff only.		
<b>Purification:</b>	The antibody was purified by affinity chromatography, and conjugated with biotin under optimal conditions. The solution is free of unconjugated biotin.		
<b>Applications:</b>	FC, ICFC		
<b>Recommended Usage:</b>	Each lot of this antibody is quality control tested by immunofluorescent staining with flow cytometric analysis as negative control. Use at concentrations comparable to those of the specific antibody of interest.		
<b>Storage &amp; Stability:</b>	The antibody solution should be stored undiluted at 4°C. <b>Do not freeze</b>		

**NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.**