

## Mouse Immunoglobulin M kappa Isotype Control Clone MM-30 LE/NA Antibody

<b>Catalog No.</b>	CSI11961 CSI11962	<b>Quantity:</b>	50 µg 500 µg
<b>Description:</b>	The isotype of MM-30 immunoglobulin is mouse IgM,k. This antibody was chosen as an isotype control after screening on a variety of resting, activated, live, and fixed mouse, rat and human tissues.		
<b>Concentration:</b>	1.0 mg/ml		
<b>Immunogen:</b>	Trinitrophenol + KLH		
<b>Isotype:</b>	Mouse IgM, k		
<b>Clone:</b>	MM-30		
<b>Formulation:</b>	Phosphate-buffered solution, pH 7.2, containing no preservative. 0.2 µm filter sterilized. Endotoxin level is < 0.1 EU/µg of the protein (< 0.01 ng/µg of the protein) as determined by the LAL test.		
<b>Purification:</b>	The LE/NA (Low Endotoxin, Azide-Free) antibody was Purified by affinity chromatography.		
<b>Applications:</b>	FC, ICFC, WB, IP, ICC, IF, FA,		
<b>Recommended Usage:</b>	Each lot of this antibody is quality control tested by immunofluorescent staining with flow cytometric analysis as negative control, and the purity is greater than 95% by SDS-PAGE. Use at concentrations comparable to those of the specific antibody of interest.		
<b>Storage &amp; Stability:</b>	The antibody solution should be stored undiluted at 4 °C. This LE/NA solution contains no preservative; handle under aseptic conditions		

**NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.**

