

Native Mouse IgG-Fc

Catalog No.	CSI11979	Quantity:	1.0 mg
Alternate Names:	Immunoglobulin G-Fc		
Description:	Native mouse IgG-Fc		
Concentration:	1.0 mg/ml		
Source:	Mouse serum - delipidated, fractionated, ion-exchange chromatography, and papain digestion		
Formulation:	Liquid in PBS, pH 7.5 with 0.01% sodium azide as preservative. Precaution: Sodium azide is a poisonous and hazardous substance which should be handled by trained staff only.		
Purity:	Single arc observed by IEP		
Storage & Stability:	Store undiluted liquid at 2-4°C. For long term storage, add glycerol to make a 50% solution and store in working aliquots at -20°C. Avoid repeated freeze-thaw cycles.		
Contaminants:	Single Band - none detected		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

