

## Recombinant Treponema pallidum p17 Protein

<b>Catalog No.</b>	CSI13231	<b>Quantity:</b>	1.0 mg
<b>Alternate Names:</b>	T. pallidum protein		
<b>Description:</b>	Purified recombinant Treponema pallidum p17 protein: Background Treponema pallidum is a species of spirochaete bacterium with subspecies that cause treponemal diseases such as syphilis, bejel, pinta and yaws. It is not seen on a Gram stained smear because the organism is too thin.		
<b>Concentration:</b>	Batch dependent - please inquire should you have specific requirements		
<b>Source:</b>	<i>E. coli</i>		
<b>Formulation:</b>	PBS with 15% glycerol and 0.1% NaN <sub>3</sub> . <b>Precaution:</b> Sodium azide is a poisonous and hazardous substance which should be handled by trained staff only.		
<b>Purity:</b>	Affinity Purified.		
<b>Species:</b>	Viral		
<b>Applications:</b>	User optimized		
<b>Storage &amp; Stability:</b>	Aliquot and store at -20 °C. <b>Avoid repeated Freeze/Thaw cycles</b>		

**NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.**