

PAK6

Recombinant Human p21 Protein-Activated Kinase 6 His-Tagged

Catalog No.	CSI13448	Quantity:	10 µg
Alternate Names:	PAK5		
Description:	Recombinant human full length PAK6, His-tagged. PAK6 is a serine/threonine kinase expressed in the testis and prostate. It associates with both the estrogen and androgen receptors. No special measures were taken to activate this kinase.		
Gene ID:	56924		
Protein Accession No:	NP_064553		
Source:	Insect cells		
Molecular Weight:	79.9 kDa		
Formulation:	Liquid in 50 mM Tris, pH 7.5 + 150 mM NaCl + 0.5 mM EDTA + 0.05% Triton X-100 + 4 mM DTT + 50% glycerol		
Storage & Stability:	Stable for 6 months in working aliquots at -80°C. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

