

KLK, Native Porcine Kallikrein

Catalog No.	CSI14903A	Quantity:	250 units
	CSI14903B		1000 units

Alternate Names: Kininogenase, Kininogenin

Description: Kallikrein is a serine protease present in blood plasma, lymph, urine, saliva, pancreatic juices, and other body fluids that catalyzes the proteolysis of plasminogen to plasmin and bradykininogen to bradykinin. In the clinical setting, Kallikreins are useful biomarkers for cancer as well as drug targets.

Source: Porcine Pancreas

Formulation: Lyophilized

Appearance: Off white to tannish powder

Biological Activity: Typically >40 U/mg powder

Specific Activity: Typically >250 U/mg Protein

Storage & Stability: Store at -20°C. 2 years.

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

