




## Western blot analysis

Mouse eye extracts and recombinant proteins (Crystallin-alpha A and B) were resolved by electrophoresis, transferred to PVDF membrane and probed with anti-Crystallin alpha A (1:1000). Proteins were visualized using a goat anti-mouse secondary antibody conjugated to HRP and a DAP detection system. This antibody is not shown cross-activity about Crystallin alpha B.



A: Crystallin-alpha A recombinant protein  
B: Crystallin-alpha B recombinant protein  
C: Mouse eye lysates

**NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.**

