

Recombinant H1N1 Influenza-A Virus (New Caledonia/20/99)

Catalog No.	CSI15838A	Quantity:	2 µg
	CSI15838B		10 µg
	CSI15838C		100 µg

Description: H1N1 is a subtype specie of Influenza A virus. H1N1 Influenza Virus has mutated into various strains such as the Spanish Flu strain, mild human flu strains, endemic pig strains, and various strains found in birds.
The Influenza A Virus is a globular particle about 100 nm in diameter, sheathed in a lipid bilayer derived from the plasma membrane of its host. Studded in the lipid bilayer are two integral membrane proteins some 500 molecules of hemagglutinin ("H") and some 100 molecules of neuraminidase ("N"). Within the lipid bilayer are 3000 molecules of matrix protein and 8 pieces of RNA. Each of the 8 RNA molecules is associated with many copies of a nucleoprotein, several molecules of the three subunits of its RNA polymerase some "non-structural" protein molecules of uncertain function.
Recombinant Full-Length H1N1 New Caledonia/20/99 is glycosylated with N-linked sugars, produced using baculovirus vectors in insect cells and its Mw is approximately 72 kDa.
The insect cells infected with A9440.1a recombinant baculovirus expressing recombinant H1N1 A/New Caledonia/20/99.
H1N1 New Caledonia shows 90% similarity to the A/PR/8/34 amino acid sequence.
The accession number is DQ508857.

Source: Baculovirus Insect Cells.

Molecular Weight: ±72 kDa

Formulation: The Recombinant H1N1 A/New Caledonia/20/99 solution 10 mM Sodium phosphate, pH 7.2 + 150 mM NaCl.

Purity: Greater than 90.0% as determined by:
Analysis by SDS-PAGE, HA1 and HA2 bands are observed using SDS-PAGE under reducing conditions.

Physical Appearance: Sterile Filtered colorless solution.

Storage & Stability: H1N1 A/New Caledonia/20/99 Recombinant should be stored at 4°C.

Applications: Western-Blot 0.1µg -1µg per strip,
ELISA 1µg/Well.

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

