

## Rabbit anti BSA IgG fraction HRP pAb

<b>Catalog No.</b>	CSI20490A CSI20490B	<b>Quantity:</b>	1.0 mg 10 mg
<b>Description:</b>	Horseradish peroxidase labeled polyclonal antibody (host rabbit). IgG fraction purified by immobilized Protein A.		
<b>Concentration:</b>	2.0 mg/ml		
<b>Specificity:</b>	Bovine		
<b>Host:</b>	Rabbit		
<b>Immunogen:</b>	Full length native protein (purified) (Cow).		
<b>Isotype:</b>	IgG		
<b>Conjugate:</b>	HRP		
<b>Formulation:</b>	Liquid in PBS, pH 7.4 + 1% BSA + 0.01% Thiomersal. Precaution: Thiomersal is a poisonous and hazardous substance which should be handled by trained staff only.		
<b>Purification:</b>	IgG fraction purified by immobilized Protein A		
<b>Cross-Reactivity:</b>	Human, Bovine		
<b>Applications:</b>	ELISA		
<b>Application Notes:</b>	The optimal concentration should be determined by the user for each specific application.		
<b>Storage &amp; Stability:</b>	When stored at -80°C, product is stable for 5 years from date of delivery. <b>Avoid repeated freeze-thaw cycles.</b>		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

