

## Rabbit Anti-Conjugated Nitrendipine pAb

<b>Catalog No.</b>	PA2069	<b>Quantity:</b>	100 µl
<b>Description:</b>	Polyclonal antisera were raised in rabbits after immunization with the Nitrendipine adsorbed on protein.		
<b>Target:</b>	Conjugated Nitrendipine		
<b>Purity:</b>	Antiserum previously pre absorbed on protein carriers, and purified.		
<b>Host:</b>	Rabbit		
<b>Immunogen:</b>	Synthetic Nitrendipine conjugated to protein carriers		
<b>Isotype:</b>	IgG		
<b>Specificity:</b>	Using a conjugate Nitrendipine-protein carrier, antibody specificity was performed with an ELISA test by competition experiments with the following compounds:		

<b>Compounds</b>	<b>Cross-reactivity ratio (a)</b>
Nitrendipine-BSA	1
Nifedipine-BSA	1/10,000

(a): Nitrendipine-BSA concentration/other conjugated compounds concentration at half displacement; BSA : Bovine Serum Albumin.

<b>Formulation:</b>	Lyophilized
<b>Purification:</b>	Antiserum previously pre-absorbed on protein carriers, and purified.
<b>Applications:</b>	ELISA
<b>Application Notes:</b>	Recommended dilutions for ELISA (1/2,000-1/5,000). Optimal dilutions should be determined by each laboratory for each application
<b>Storage &amp; Stability:</b>	Lyophilized vial must be stored at 4°C in a dry area. After reconstitution with 50 µl of distilled water and 50 µl of glycerol, the aliquot can be stored at -20°C, and is stable at least 2 years.
<b>Research Areas:</b>	Neurobiology

**NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.**

