

Rabbit Anti-Morganella Morganii pAb

Catalog No.	PA2072	Quantity:	100 µl
Target:	<i>Morganella morganii</i>		
Purity:	Ammonium Sulphate purified		
Host:	Rabbit		
Immunogen:	<i>Morganella morganii</i> total proteins		
Isotype:	IgG		
Specificity:	Antibody specificity was performed with an ELISA test by competition experiments with the following compounds :		

Compounds	Cross-reactivity ratio (a)
<i>Morganella morganii</i>	1
<i>Hafnia alvei</i>	1/50
<i>Pseudomonas aeruginosa</i>	1/50
<i>Pseudomonas putida</i>	1/50
<i>Proteus mirabilis</i>	1/50
<i>Citrobacter koseri</i>	1/50
<i>Klebsiella pneumoniae</i>	1/50

(a): *Morganella morganii* antigens concentration / Other bacteria antigens concentration at half displacement

Formulation:	Lyophilized
Purification:	Antiserum previously purified
Applications:	ELISA
Application Notes:	Recommended dilutions for Western Blot (1/1,000-1/2,000)
Storage & Stability:	Lyophilized vial must be stored at 4°C in a dry area. After reconstitution with 50 µl of distilled water and 50 µl of glycerol, the aliquot can be stored at -20°C, and is stable at least 2 years.
Research Areas:	Bacteriology, Infectious diseases, Chronic diseases
Notes:	Principal name <i>Morganella morganii</i>

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

