


## LGALS9 Human Galectin-9 ELISA Kit

<b>Catalog No.</b>	CDK340A CDK340B	<b>Quantity:</b>	1 x 96 tests 2 x 96 tests
<b>Alternate Names:</b>	Gal-9, Ecalectin, Tumor antigen HOM-HD-21, HUAT		
<b>Description:</b>	The Cell Sciences™ Galectin-9 ELISA kit is a solid phase sandwich ELISA for the in-vitro qualitative and quantitative determination of Galectin-9 in supernatants, buffered solutions or serum and plasma samples. This assay will recognize both native and recombinant human Galectin-9.		
<b>UniProt ID:</b>	O00182		
<b>Range:</b>	500 - 15.6 pg/ml		
<b>Sensitivity:</b>	8 pg/ml		
<b>Sample Size:</b>	50 µl or 100 µl		
<b>Incubation Time:</b>	3 hours 45 minutes		
<b>Precision:</b>	Intra-assay CV 3.2% Inter-assay CV 5.7%		
<b>Recovery:</b>	Serum, Media Range 79 - 120% Mean 99%		
<b>Contents:</b>	ELISA Kit includes pre-coated strip plate, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents and a detailed procedure.		
<b>Storage &amp; Stability:</b>	Store entire kit at 2-8°C.		
<b>Cross-Reactivity:</b>	No cross reactivity with GM-CSF, MICB, IL-34, IFNγ, IL-29, Galectin-1, Galectin-3, Galectin-8, Galectin-10		



NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.