


IL15

Human IL-15 ELISA Kit

Catalog No.	CDK341A CDK341B	Quantity:	1 x 96 tests 2 x 96 tests
Alternate Names:	Interleukin-15,		
Description:	The Cell Sciences™ IL-15 ELISA kit is a solid phase sandwich ELISA for the in-vitro qualitative and quantitative determination of IL-15 in supernatants, buffered solutions or serum and plasma samples. This assay will recognize both native and recombinant human IL-15.		
UniProt ID:	P40933		
Range:	1000 - 31.25 pg/ml		
Sensitivity:	18.3 pg/ml		
Sample Size:	50 µl and 100 µl		
Incubation Time:	3 hours 45 minutes		
Precision:	Intra-assay CV 3.6% Inter-assay CV 3.8%		
Recovery:	Serum Range 64 - 86% Mean 77%		
Contents:	ELISA Kit includes pre-coated strip plate, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents and a detailed procedure.		
Storage & Stability:	Store entire kit at 2-8°C.		
Cross-Reactivity:	No cross reactivity with IL-2, IL-4, IL-7, IL-9, IL-10, IL-13, IL-21 and IFN-g.		



NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

