

Csf2

Mouse GM-CSF ELISA Kit

Catalog No. CDK344A **Quantity:** 1 x 96 tests
CDK344B 2 x 96 tests

Alternate Names: Granulocyte macrophage colony stimulating factor

Description: The Cell Sciences™ GM-CSF ELISA kit is a solid phase sandwich ELISA for the in-vitro qualitative and quantitative determination of murine GM-CSF in supernatants, buffered solutions or serum and plasma samples. This assay will recognize both native and recombinant murine GM-CSF.

UniProt ID: P01587

Range: 500 - 15.6 pg/ml

Sensitivity: < 9 pg/ml

Sample Size: 100 µl


Incubation Time: 3 hours 45 minutes

Precision: Intra-assay CV 4.6%
Inter-assay CV 5.2%

Contents: ELISA Kit includes pre-coated strip plate, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents and a detailed procedure.

Storage & Stability: Store entire kit at 2-8°C.

Cross-Reactivity: No cross reactivity with other murine cytokines.



NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

