


## IL6

### Mouse IL-6 ELISA Kit

<b>Catalog No.</b>	CDK347A CDK347B	<b>Quantity:</b>	1 x 96 tests 2 x 96 tests
<b>Alternate Names:</b>	Interleukin-6, B-cell hybridoma growth factor, Interleukin HP-1		
<b>Description:</b>	The Cell Sciences™ IL-6 ELISA kit is a solid phase sandwich ELISA for the in-vitro qualitative and quantitative determination of murine IL-6 in supernatants, buffered solutions or serum and plasma samples.		
<b>UniProt ID:</b>	P08505		
<b>Range:</b>	500 - 15.6 pg/ml		
<b>Sensitivity:</b>	< 10 pg/ml		
<b>Specificity:</b>	This assay will recognize both native and recombinant murine IL-6.		
<b>Sample Size:</b>	100 µl		
<b>Incubation Time:</b>	3 hours 45 minutes		
<b>Precision:</b>	Intra-assay CV 4.7% Inter-assay CV 4.4%		
<b>Contents:</b>	ELISA Kit includes pre-coated strip plate, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents and a detailed procedure.		
<b>Storage &amp; Stability:</b>	Store entire kit at 2-8°C.		



NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

