


Ifng

## Rat IFN-gamma ELISA Kit

<b>Catalog No.</b>	CDK349A CDK349B	<b>Quantity:</b>	1 x 96 tests 2 x 96 tests
<b>Alternate Names:</b>	Interferon gamma, IFN- $\gamma$		
<b>Description:</b>	The Cell Sciences™ IFN-g ELISA kit is a solid phase sandwich ELISA for the in-vitro qualitative and quantitative determination of rat IFN-g in supernatants, buffered solutions or serum and plasma samples.		
<b>UniProt ID:</b>	P01581		
<b>Range:</b>	1000 - 31.25 pg/ml		
<b>Sensitivity:</b>	< 10 pg/ml		
<b>Specificity:</b>	Recognizes both native and recombinant rat IFN-gamma		
<b>Sample Size:</b>	100 $\mu$ l		
<b>Incubation Time:</b>	3 hours 45 minutes		
<b>Precision:</b>	Intra-assay CV 3.2% Inter-assay CV 4.9%		
<b>Contents:</b>	ELISA Kit includes pre-coated strip plate, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents and a detailed procedure.		
<b>Storage &amp; Stability:</b>	Store entire kit at 2-8°C.		
<b>Cross-Reactivity:</b>	No cross reactivity with other rat cytokines.		



NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

