


IL4

Rat IL-4 ELISA Kit

Catalog No.	CDK350A CDK350B	Quantity:	1 x 96 tests 2 x 96 tests
Alternate Names:	Interleukin-4, B-cell growth factor 1, B-cell stimulatory factor 1, BSF-1, Lymphocyte stimulatory factor 1		
Description:	The Cell Sciences™ IL-4 kit is a solid phase sandwich ELISA for the in-vitro qualitative and quantitative determination of rat IL-4 in supernatants, buffered solutions or serum and plasma samples.		
UniProt ID:	P20096		
Range:	1000 - 31.25 pg/ml		
Sensitivity:	< 15 pg/ml		
Specificity:	Recognizes both native and recombinant rat IL-4		
Sample Size:	100 µl		
Incubation Time:	4 hours 45 minutes		
Precision:	Intra-assay CV 2.1% Inter-assay CV 4.9%		
Contents:	ELISA Kit includes pre-coated strip plate, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents and a detailed procedure.		
Storage & Stability:	Store entire kit at 2-8°C.		
Cross-Reactivity:	No cross reactivity with other rat cytokines.		



NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

