


PTPRC

Mouse Anti-Human CD45 (Clone B-A11) mAb

Catalog No.	CDM188	Quantity:	200 tests
Alternate Names:	LCA; LY5; B220; CD45; L-CA; T200; CD45R; GP180		
Description:	The mouse monoclonal antibody recognizes human CD45, a member of the protein tyrosine phosphatase (PTP) family. PTPs are known to be signaling molecules that regulate a variety of cellular processes including cell growth, differentiation, mitosis, and oncogenic transformation. CD45 is a receptor type PTP, with an extracellular domain, a single transmembrane segment and two tandem intracytoplasmic catalytic domains. CD45 is an essential regulator of T- and B-cell antigen receptor signaling. It functions through either direct interaction with components of the antigen receptor complexes, or by activating various Src family kinases required for the antigen receptor signaling. CD45 also suppresses JAK kinases, and thus functions as a regulator of cytokine receptor signaling.		
Gene ID:	5788		
Concentration:	100 tests/ml		
Conjugate:	Unconjugated		
Specificity:	Recognizes the Leucocyte Common Antigen (LCA) with different isoforms of 180-190-205-220 kDa.		
Host:	Mouse		
Isotype:	IgG1		
Immunogen:	Thymus cells and Jurkat cell line.		
Clone:	B-A11		
Formulation:	Liquid in PBS + 1 % BSA + 0.09 % sodium azide. Precaution: Sodium azide is a poisonous and hazardous substance which should be handled by trained staff only.		
Applications:	Flow Cytometry, Immunohistochemistry		
Application Notes:	Use 10 µl to label 10 ⁶ cells or 100 µl of whole blood. The optimal concentration should be determined by the user for each specific application.		
Storage & Stability:	Stable at 2-8°C for 1 year. For longer storage, freeze aliquots at -20 to -80°C. Avoid repeated freeze-thaw cycles.		

A typical staining pattern with the B-A11 monoclonal antibody of lymphocytes



NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences[®]
65 Parker Street
Unit 11
Newburyport, MA 01950

Toll Free: 888-769-1246
Phone: 978-572-1070
Fax: 978-992-0298

E-mail: info@cellsciences.com
Website: www.cellsciences.com