


## PTPRC

### Mouse Anti-Human CD45RA (Clone B-C15) mAb

<b>Catalog No.</b>	CDM192	<b>Quantity:</b>	200 tests
<b>Alternate Names:</b>	Receptor-type tyrosine-protein phosphatase C, Leukocyte common antigen, L-CA, T200, CD45		
<b>Description:</b>	The CD45RA antigen is a member of the CD45 antigen family that also includes the CD45, CD45RB, and CD45RO antigens. The CD45RA antigen is present on approximately 50% of CD4+ T lymphocytes, approximately 75% of CD8+ T lymphocytes, and on essentially all B lymphocytes and natural killer lymphocytes. <sup>19</sup> The suppressor inducer T-lymphocyte subset expresses the phenotype CD4+CD45RA+. <sup>19</sup> The CD45RA antigen is expressed on naive T lymphocytes; antigen density decreases upon in vitro activation. <sup>20</sup> A selective loss of the CD4+CD45RA+ subset during active multiple sclerosis has been demonstrated.		
<b>UniProt ID:</b>	P08575		
<b>Gene ID:</b>	5788		
<b>Hybridoma:</b>	Myeloma X63/AG.8653 x BALB/c spleen cells		
<b>Specificity:</b>	Recognizes the restricted Leucocyte Antigen (LCA) with different isoforms of 205-220 kDa		
<b>Host:</b>	Mouse		
<b>Immunogen:</b>	Reh6 cell line		
<b>Isotype:</b>	IgG1		
<b>Clone:</b>	B-C15		
<b>Concentration:</b>	200 tests / 2 ml		
<b>Formulation:</b>	PBS, 1% BSA, 0.1% sodium azide. Precaution: Sodium azide is a poisonous and hazardous substance which should be handled by trained staff only.		
<b>Purification:</b>	Ion exchange chromatography		
<b>Applications:</b>	Flow Cytometry		
<b>Application Notes:</b>	Use 10 µl to label 10 <sup>6</sup> cells or 100 µl of whole blood.		
<b>Storage &amp; Stability:</b>	Stable at 2-8°C for 12 months. For longer storage, freeze aliquots. <b>Avoid repeated freeze-thaw cycles.</b>		

A typical staining pattern with the B-C15 monoclonal antibody of lymphocytes.



NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



**Cell Sciences**<sup>®</sup>  
65 Parker Street  
Unit 11  
Newburyport, MA 01950

Toll Free: 888-769-1246  
Phone: 978-572-1070  
Fax: 978-992-0298

E-mail: [info@cellsciences.com](mailto:info@cellsciences.com)  
Website: [www.cellsciences.com](http://www.cellsciences.com)