


TEK

Mouse Anti-Human Tyrosine Endothelial Kinase/TIE-2 Clone B-B48 mAb

Catalog No.	CDM229	Quantity:	200 µg
Alternate Names:	CD202B		
Description:	Mouse Anti-Human TEK/TIE-2, Clone B-B48, Monoclonal Antibody		
Concentration:	200 µg/200 µl		
Specificity:	TEK/TIE-2		
Host:	Mouse		
Immunogen:	Eahy 926 cell line		
Isotype:	Mouse IgG1		
Clone:	B-B48		
Formulation:	Liquid sterile filtered in PBS, treated to remove endotoxins. Carrier and preservative free.		
Purification:	Ion exchange chromatography		
Applications:	Suitable for flow cytometry Inhibits Tek-L / Tek binding		
Storage & Stability:	Stable for 1 year at 2-8°C. For longer storage, freeze aliquots at -20°C. Avoid repeated freeze-thaw cycles.		

A typical staining pattern with the B-B48 monoclonal antibody of Eahy 926 cell line.



NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.