

Mouse IgG2b Isotype Control (Clone B-F33) mAb, Azide Free

Catalog No.	CDM513A CDM513B	Quantity:	200 µg 500 µg
Description:	The mouse monoclonal antibody is for use as an <i>in vivo</i> isotype control.		
Hybridoma:	Myeloma X63/AG.8653 x Balb/c spleen cells		
Concentration:	1 mg/ml		
Conjugate:	Unconjugated		
Specificity:	Not specific for any known human antigens.		
Host:	Mouse		
Isotype:	IgG2b λ		
Immunogen:	Chemical molecule		
Clone:	B-F33		
Formulation:	Sterile-filtered PBS, treated to remove endotoxins. Carrier and preservative free.		
Applications:	Flow Cytometry, IHC, <i>in vivo</i> control for functional assays.		
Storage & Stability:	Stable at 2-8°C for 1 year. DO NOT FREEZE.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

