


Mouse Anti-Goat IgG Fc Clone SB115d HRP mAb

Catalog No.	CMG002	Quantity:	1.0 mL
Description:	Mouse Anti-Goat IgG Fc Clone SB115d HRP mAb		
Specificity:	Goat IgG Fc		
Immunogen:	Pooled goat IgG		
Isotype:	Mouse (BALB/c) IgG _{2a} K		
Clone:	SB115d		
Conjugate:	HRP (Horseradish Peroxidase)		
Formulation:	Liquid volume of 1.0 mL in a stock solution in 50% glycerol/50% PBS, pH 7.4. No preservative added.		
Isotype Control	Mouse IgG2a-HRP		
Cross-Reactivity:	Minimal reactivity to goat IgG Fab and sheep, rabbit, human, monkey, rat, hamster, canine, and swine immunoglobulins		
Applications:	Applications for relevant formats of this clone include; ELISA, and FLISA.		
Application Notes:	Working dilutions; ELISA conjugate 1:5,000 – 1:20,000. Since applications vary, you should determine the optimum working dilution for the product that is appropriate for your specific need.		
Storage & Stability:	Store at 2-8°C or long-term at -20°C.		

ELISA plate was coated with purified mouse IgG Fc and IgG Fab and human, bovine, and horse serum. Immunoglobulins and sera were detected with serially diluted Goat Anti-Mouse IgG Fc, Human/Bovine/Horse SP ads-UNLB followed by Mouse Anti-Goat IgG Fc-HRP.



NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

