


Mouse Anti-Human Kappa Clone SB81a HRP mAb

| | | | |
|---------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------|--------|
| Catalog No. | CMH035 | Quantity: | 1.0 mL |
| Description: | Mouse Anti-Human Kappa Clone SB81a HRP mAb | | |
| Specificity: | Human Kappa | | |
| Isotype: | Mouse (BALB/c) IgG1κ | | |
| Clone: | SB81a | | |
| Conjugate: | HRP (Horseradish Peroxidase) | | |
| Formulation: | Liquid volume of 1.0 mL in a stock solution of 50% glycerol/50% PBS, pH 7.4. No preservative added. | | |
| Isotype Control | Mouse IgG ₁ -HRP | | |
| Applications: | Applications for relevant formats of this clone - ELISA. | | |
| Application Notes: | Working dilutions for ELISA - 1:4,000 - 1:8,000. Since applications vary, you should determine the optimum working dilution for the product that is appropriate for your specific need. | | |
| Storage & Stability: | Store at 2-8°C or long-term at -20°C. | | |

ELISA plate was coated with purified human IgGκ, IgMκ, IgAκ, IgGλ, IgMλ, and IgAλ. Immunoglobulins were detected with serially diluted Mouse Anti-Human Kappa-HRP.



NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

