

Ifng

Rat Anti-Mouse Interferon gamma (Clone F3) mAb

Catalog No.	CMI035	Quantity:	0.5 mg
Alternate Names:	IFN-gamma		
Description:	Rat Anti-Mouse Interferon gamma (Clone F3) is a neutralizing monoclonal antibody.		
UniProt ID:	P01580		
Gene ID:	15978		
Source:	Produced <i>in vitro</i> using serum free medium.		
Specificity:	Binds with high efficacy to Native and Recombinant Mouse IFN gamma. Neutralizing <i>In vivo</i> and <i>in vitro</i> .		
Isotype:	Rat IgG _{2a}		
Clone:	F3		
Formulation:	Lyophilized from sterile-filtered PBS containing 125 mM Trehalose		
Purification:	Ion Exchange Chromatography		
Reconstitution:	Centrifuge vial briefly before opening. Reconstitute with 0.5 ml sterile distilled water, containing 0.02% sodium azide to prevent bacterial growth (recommended). Precaution: Sodium azide is a poisonous and hazardous substance which should be handled by trained staff only.		
Applications:	<i>In vivo</i> and <i>in vitro</i> neutralization ELISA, ELISPOT assays		
Storage & Stability:	Lyophilized product is stable for at least one year at 2-8°C. After reconstitution, the contents can be safely stored at 2-8°C for one month or for one year in working aliquots at -20°C to -80°C. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences[®]
65 Parker Street
Unit 11
Newburyport, MA 01950

Toll Free: 888-769-1246
Phone: 978-572-1070
Fax: 978-992-0298

E-mail: info@cellsciences.com
Website: www.cellsciences.com