

IL13

Mouse Anti-Human IL-13 (Clone B-B13) Cy5-labeled mAb

Catalog No.	CMI259	Quantity:	50 µg
Alternate Names:	Interleukin-13		
Description:	The mouse anti-human IL-13 monoclonal antibody (clone B-B13) can be used in ELISA, ELISPOT and for <i>in vitro</i> neutralization assays. Demonstrated to block IL-13 stimulated TARC production in human lung epithelial A-549 cells.		
UniProt ID:	P35225		
Gene ID:	3596		
Specificity:	Binds with high affinity both native and recombinant human IL-13.		
Source:	Produced <i>in vitro</i> using serum free medium.		
Isotype:	Mouse IgG1		
Clone:	B-B13		
Conjugate:	Cy5 Bis-N-hydroxysuccinimide ester		
Quantitation:	$E^{0.1\%}_{280\text{ nm}} = 1.4$ for a 1 mg/ml solution.		
Formulation:	Lyophilized from sterile-filtered PBS containing 1% BSA and 125 mM Trehalose		
Purification:	Ion exchange chromatography		
Reconstitution:	Centrifuge vial briefly before opening. Reconstitute with 0.5 ml sterile distilled water, containing 0.02% sodium azide to prevent bacterial growth (recommended). Precaution: Sodium azide is a poisonous and hazardous substance which should be handled by trained staff only.		
Applications:	Flow cytometry		
Storage & Stability:	Lyophilized product is stable for at least one year at 2-8°C. After reconstitution, the contents can be safely stored at 2-8°C for one month or for one year in working aliquots at -20°C to -80°C. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences®
65 Parker Street
Unit 11
Newburyport, MA 01950

Toll Free: 888-769-1246
Phone: 978-572-1070
Fax: 978-992-0298

E-mail: info@cellsciences.com
Website: www.cellsciences.com