

## Rat Anti-Mouse IgE Clone 23G3 HRP mAb

<b>Catalog No.</b>	CMM016	<b>Quantity:</b>	1.0 mL
<b>Description:</b>	Rat Anti-Mouse IgE Clone 23G3 HRP mAb		
<b>Specificity:</b>	Mouse IgE		
<b>Immunogen:</b>	Monoclonal anti-DNP mouse IgE		
<b>Isotype:</b>	Rat IgG <sub>1</sub> κ		
<b>Clone:</b>	23G3		
<b>Conjugate:</b>	HRP (Horseradish Peroxidase)		
<b>Formulation:</b>	Liquid volume of 1.0 mL in a stock solution of 50% glycerol/50% PBS, pH 7.4. No preservative added		
<b>Isotype Control</b>	Rat IgG <sub>1</sub> -HRP		
<b>Applications:</b>	Applications for relevant formats of this clone include ELISA, Immunohistochemistry, Western Blot and ELISPOT.		
<b>Application Notes:</b>	Working dilutions for ELISA 1:1,000 - 1:2,000. Since applications vary, you should determine the optimum working dilution for the product that is appropriate for your specific need.		
<b>Storage &amp; Stability:</b>	Store at 2-8°C or long-term at -20°C.		

**NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.**

