


Rat Anti-Mouse IgG_{2a} Clone SB84a FITC mAb

Catalog No.	CMM024	Quantity:	500 µg
Description:	Rat Anti-Mouse IgG _{2a} Clone SB84a FITC mAb		
Concentration:	500 µg/mL		
Specificity:	Mouse IgG _{2a}		
Immunogen:	Mouse IgG _{2a} hybridoma		
Isotype:	Rat (Sprague Dawley) IgG _{1κ}		
Clone:	SB84a		
Conjugate:	FITC (Fluorescein)		
Formulation:	Liquid volume in 1.0 mL in a buffer of PBS/NaN ₃		
Isotype Control	Rat IgG ₁ -FITC		
Applications:	Applications for relevant formats of this clone include ELISA, Flow Cytometry, Immunohistochemistry, Immunocytochemistry, Western Blot, and Multiplex.		
Application Notes:	Working dilutions; FLISA 1:200 - 1:400. Since applications vary, you should determine the optimum working dilution for the product that is appropriate for your specific need.		
Storage & Stability:	Store at 2-8°C. DO NOT FREEZE		

FLISA plate was coated with purified mouse IgG1, IgG2a, IgG2b, IgG3, IgM, and IgA. Immunoglobulins were detected with serially diluted Rat Anti-Mouse IgG2a-FITC.



NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.