


## Rat Anti-Mouse Kappa Clone 187.1 HRP mAb

<b>Catalog No.</b>	CMM030	<b>Quantity:</b>	1.0 mL
<b>Description:</b>	Rat Anti-Mouse Kappa Clone 187.1 HRP mAb		
<b>Specificity:</b>	Mouse Kappa		
<b>Immunogen:</b>	MPC-11 immunoglobulin		
<b>Isotype:</b>	Rat (Sprague Dawley) IgG <sub>1</sub> κ		
<b>Clone:</b>	187.1		
<b>Conjugate:</b>	HRP (Horseradish Peroxidase)		
<b>Formulation:</b>	Liquid volume of 1.0 mL supplied in a stock solution in 50% glycerol/50% PBS, pH 7.4. No preservative added		
<b>Isotype Control</b>	Rat IgG <sub>1</sub> -HRP		
<b>Applications:</b>	Applications for relevant formats of this clone include ELISA, Immunohistochemistry, and Western Blot.		
<b>Application Notes:</b>	Recommended working dilution for ELISA is 1:4,000-1:8,000. Since applications vary, you should determine the optimum working dilution for the product that is appropriate for your specific need.		
<b>Storage &amp; Stability:</b>	Store at 2-8°C or long-term at -20°C.		

ELISA plate was coated with purified mouse IgGκ, IgMκ, IgAκ, IgGλ, IgMλ, and IgAλ. Immunoglobulins were detected with serially diluted Rat Anti-Mouse Kappa-HRP.



NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

